

ANEXdb

Animal Expression Database

USER'S MANUAL



Revision Sheet

Release No.	Date	Author	Revision Description
Rev. 0	12/15/08	Sushain Pandit	User's Manual First-cut
Rev. 1	12/25/08	Sushain Pandit	Added Section 3
Rev. 2	01/05/09	Sushain Pandit	Added Sections 4 & 5

USER'S MANUAL

TABLE OF CONTENTS

	<u>Page #</u>
1.0 GENERAL INFORMATION	1-1
1.1 System Overview	1-1
1.2 References	1-1
1.3 Authorized Use Permission	1-1
1.4 Points of Contact	1-1
1.5 Acronyms and Abbreviations	1-1
2.0 GETTING STARTED	2-1
2.1 Logging On.....	2-1
2.2 Dashboard.....	2-2
2.3 Exiting the System	2-3
3.0 USING THE SYSTEM.....	3-1
3.1 Submitting an experiment	3-1
3.2 Uploading associated files.....	3-2
3.3 Hybridizations page.....	3-3
3.4 Review page	3-4
4.0 QUERYING	4-1
4.1 Query Capabilities.....	4-1
5.0 IMPORTANT INSTRUCTIONS AND FAQs.....	5-1
5.1 Submission and Finalization process	5-1
5.2 Regarding SOFT files and zip archives.....	5-2
5.3 Submitting the experiment to GEO.....	5-2
5.4 Updating the SOFT after submission (if required).	5-3

1.0 GENERAL INFORMATION

1.0 GENERAL INFORMATION

1.1 System Overview

ANEXdb allows scientists to store and analyze gene expression data in the form of microarray experiments. Currently [Affymetrix](#) is the only platform supported. The ANEXdb database schema is [available here](#).

- The underlying functions performed by the system are storage and query-based retrieval.
- The system is *broadly* based on a web-based, 3-tier architecture.
- User access mode is essentially restricted via the web-based client.
- System status is *Operational*.

1.2 References

This document was prepared as is by the ANEXdb core development team. Refer www.anexdb.org.

1.3 Authorized Use Permission

This system is protected from unauthorized usage and/or making unauthorized copies of data, software pieces, etc. For further information on this, please get in touch with the point of contact(s) (Refer 1.4).

1.4 Points of Contact

Please drop a note to simius@iastate.edu for queries or feedback.

1.5 Acronyms and Abbreviations

ANEXdb	Animal Expression Database
NSF	National Science Foundation
USDA	United States Department of Agriculture

2.0 GETTING STARTED

2.0 GETTING STARTED

2.1 Logging On / Registering

Depending on whether the system is being used for the first time, you may need to register by clicking 'Register' at the top-right part on the landing page at www.anexdb.org .

ANEXdb - Register - Mozilla Firefox

File Edit View History Bookmarks Tools Help

http://www.anexdb.org/register.php

Getting Started Latest Headlines svn cat ANEXdb - Animal Expr...

ANEXdb - Animal Expression ... ANEXdb - Animal Expression ... ANEXdb - Animal Expression ... ANEXdb - Animal Expression ... ANEXdb - Animal Expression ... ANEXdb - Register

ANEXdb

ANEXdb: Animal Expression Database

Home Search Help Register | Login

ANEXdb > Register

Register

*Email address:

*Re-enter email address:

First name:

Last name:

Job title:

Head of lab:

Lab name:

Institution:

Street address:

City:

State:

Country:

Zip Code:

Done

start anexdb.org - default... anexdb.org - default... Java - anexdb/src/ed... user_manual - Micro... ANEXdb - Register - ... 12:34 PM

Filling in all the fields (besides mandatory) is recommended so that the admin can track the submissions. When you register, it is only required that you enter an email address and password, but the more of the other details you enter here, the less you will have to enter later when submitting sample information.

2.2 Dashboard

After you register, login using your email address and password. After logging into the system, the user will be presented with a dashboard, where he can create new experiments.

The screenshot shows the ANEXdb Dashboard interface. At the top, there is a navigation bar with links for Home, Dashboard, Search, and Help. The user is logged in as 'sushain.pandit@gmail.com' and has a 'Logout' link. The main content area is divided into several sections:

- Your Groups:** A table with columns 'Group Name' and 'Created'. Below the table, it states 'You are not currently a member of any groups.'
- Your Account Info (edit):** A form displaying user details:

Email	sushain.pandit@gmail.com
FirstName	Sushain
LastName	Pandit
JobTitle	Graduate Research Assistant
HeadOfLab	Dr.Honavar
Lab	Artificial Intelligence Research Lab
Institution	Iowa State University
Street	
State	
City	
Country	
Zip	
Phone	515-451-0757
- Your Submissions:** A table with columns 'Experiment', 'Date Created', and 'Status'. Below the table, it states 'You have no pending experiment submissions.'
- Your Experiments:** A table with columns 'Experiment', 'Date Created', and 'visibility'. It contains the following data:

Experiment	Date Created	visibility
yet_another_test	2009-01-02 00:00:00	Public
another_test_2	2009-01-02 00:00:00	Public
testing_multiple_hybrids	2009-01-02 00:00:00	Public
testing_new_jar	2009-01-02 00:00:00	Public

At the bottom of the dashboard, there is a footer with the text 'Questions or comments: simius@iastate.edu'.

2.2.1 Start New Submission

This may be used to start the new submission. Click on the "Start New Submission" link to begin entering a new experiment.

2.2.1 Your Submissions / Experiments

These boxes enlist all the submissions done under one particular user-id along with other useful information on the current status of submissions.

2.2.1 Account Info / Groups

These enlist your profile information and whether you've any group privileges as determined by the ANEXdb admin.

2.3 Exiting the System

Make sure you logout properly by clicking 'logout' on the top right. Any incomplete submissions can be located within 'Your Submissions' box when you log into the system the next time.

3.0 USING THE SYSTEM

3.0 USING THE SYSTEM

3.1 Submitting an experiment

Clicking ‘Start New Submission’ on your dashboard starts a new experiment.

Experiment
ANEXdb follows the guidelines for a GEO SOFT submission. For a full description go here.

*Experiment Title:

*Experiment Description:

*Visibility: Group (optional):

GEO Accession Number (GSE or GDS):

Experiment Design

Factor 1: Number of Levels:

Factors should include every factor relevant to the study:

- Experimental conditions such as time points, animal lines, administration of a drug and dose
- Animal specific items such as animal numbers
- Replicate information for either biological or technical replicates

Following is a general guideline for experiment submissions.

Title: The experiment title is the same as what you would see in GEO. An example would be "Expression data from non-infected and Salmonella Typhimurium infected mesenteric lymph nodes".

Description: This is a longer text description of the experiment and its design. You should probably mention the factors and methods used in the experiment here.

Visibility: Here you will set whether you would like the experiment to be available to the public. Public visibility is the default setting and allows anyone to view your experiment and data. Private keeps the experiment hidden so that only you can view it. If you are a member of a Group you will also be able to choose group visibility and that group so that only you and other members of your group can view the experiment.


GEO accession number (optional): If your experiment has already been submitted and reviewed by NCBI GEO, you can enter the GDS or GSE accession number here.

Experiment Design: Enter the Factors and Levels of your experiment here. For example, if your experiment is a time course experiment with four time steps, you would enter "Time" in the Factor box, and "4" in the Number of Levels box. If you have only this factor, you would then click "Refresh Factors" to enable 4 text boxes so that you could enter your levels (e.g. 4 hours, 8 hours, 16 hours, 24 hours).

Publications (optional): If a journal article related to this experiment has already been published, you can enter the article details here.

Links (optional): If you have a website with supplementary material about this experiment, you can enter that link here.

3.2 Uploading associated files



Once the plugin loads, click "Browse..." and select the .CEL files and .CHP files from your experiment. You should probably have some way of identifying by the filename which chip it belongs to so that you can associate it on the next page. When you are done selecting the files, be sure to click the "Upload" button. Once the upload is done, move on to the next page.

3.3 Hybridizations page

Hybridization
ANEXdb follows the guidelines for a GEO SOFT submission. For a full description go here.

Select Array Design: Choose an Array Design ?

GEO Accession Number (GSM): ?

Select CEL File: Choose a CEL File ?

Select CHP File: Choose a CHP File ?

Factors
?

Sample

Select Organism: Select an Organism ?

Organism Part: Part 1 ?

Sample Molecule: total RNA ?

Labeling Compound: ?

Sample Description: ?

Treatment Protocol: ?

Following is a general guideline for entering the data on this page.

Array Design: This drop-down menu lists all of the Affymetrix chips currently loaded. If your experiment used a chip that is not listed, please contact an administrator using the email address at the bottom of this page.

GEO Accession Number: If your experiment has already been submitted and reviewed by NCBI GEO, you can enter the GSM accession number associated with this particular chip here.

CEL File: This drop-down will contain a list of the CEL files you uploaded on the previous page. Select the one for which you would like to add details.

CHP file: This drop-down will contain a list of the CHP files you uploaded on the previous page. Select the one for which you would like to add details.

Factors: Select the levels of the factors related to the chip you selected. If you selected pig_4hours.cel, select "4 hours" in the "Time" drop-down menu here.

Organism: Select the organism from which your sample came.

Organism part: Enter the type of tissue from which the sample came. Examples are liver, lung, MLN.

Sample molecule: The molecule type extracted from the sample.

Labeling compound: This is usually biotin for Affymetrix experiments. In dual-color experiments this is often Cy3 or Cy5.

Sample description: A short description of this particular sample. An example would be "mesenteric lymph nodes, non-infected pig, rep 1".

Treatment protocol: Describe the sample treatment. You can find an example of this at [this link](#).

Extract, Scan, Label, Hybrid Protocols: These drop-down menus allow you to quickly select a complete description of your protocol. If you do not find your particular protocol in the list, please contact us.

When you are done with that page, you have several options: -

- a) Cancel - This will not submit your current hybridization and will bring you back to the main page.
- b) Save - This will save the hybridization you are working on, but will keep you on that page.
- c) Add another hybridization - Saves your current hybridization and lets you add another .
- d) Save & Continue - If this is your last hybridization, click this to save and move on to the last review step.
- e) Continue without Saving - If this is your last hybridization, click this to move on to the last review step without saving this hybridization.


3.4 Review page

This last page allows you to look at what you've submitted and make sure it's correct. When you are ready to submit, click the Complete Submission button and it will be approved by an administrator before being submitted to the public database. You should now see your experiment in the Your Submission table in "Pending" status.

4.0 QUERYING

4.0 QUERYING

The system can be queried using varied criteria on the search page - <http://www.anexdb.org/search.php>



4.1 Query Capabilities

4.1.1 Search by Experiment

This option can be used to search using any of the following:

- Experiment Accession number (ANX.species.#)
- Experiment title keywords
- GEO Data Set or Series (GDS/GSE)

4.1.2 Search by Consensus ID

This option can be used to search using:

- Enter Cluster number followed by Contig number, i.e. ITC.1.1 - GEO Data Set or Series (GDS/GSE)

4.1.3 Search by Hybridization

This option can be used to search using any of the following:

- Hybridization accession number (ANX.species.#.#)
- GEO Sample Number (GSM)
- Sample description
- Sample tissue

4.1.4 Search by Gene or Sequence ID

This option can be used to search using any of the following:

- Gene Name
- Affymetrix Id
- NCBI RefSeq ID

4.1.5 Search by Affymetrix Probe ID

This option can be used to search using affymetrix probe id. For instance - Ssc.10082.1.A1_at.

4.1.6 Browse by Organism

This option can be used to search using any of the organism types using the drop-down menu given.

5.0 IMPORTANT INSTRUCTIONS AND FAQs

5.0 IMPORTANT INSTRUCTIONS AND FAQs

This section contains some important observations that may come handy to the users.

5.1 Submission and Finalization process

- Once an experiment is submitted, it is temporarily stored in an intermediate database and remains there until the admin *finalizes* it.
- The finalization process is essentially controlled by the system admin and it is his responsibility to ensure that the experiment submission is proper and in accordance to guidelines, however, to make his job easier, it is recommended that submissions be made as per the above specified guidelines.
- Once the experiment gets finalized, it is queued and gets processed within a pre-determined time-slot by a server-based *cron* process. This process runs 12 am – 7 am due to server load constraints and thus, it is not recommended to browse / click / review your experiment until a day after it is finalized and appears in the Search Results. This will ensure that the data has been processed and the soft files will be properly generated.

Failure to conform with the above guideline will result in incorrectly populated soft file references.

- Once the experiment is finalized (can be viewed in the Search Results) and the above mentioned time-period has elapsed, one may go ahead and click on the corresponding experiment on the search page.

The screenshot displays the ANEXdb Animal Expression Database search results page. The page title is "ANEXdb: Animal Expression Database" and the user is logged in as "sushain.pandit@gmail.com". The search results are displayed in a table with the following columns: Experiment, Organism, Contact, and Date Created.

Experiment	Organism	Contact	Date Created
yet_another_test	Bos taurus	Sushain Pandit	2009-01-02 00:00:00
another_test_2	Ovis aries	Sushain Pandit	2009-01-02 00:00:00
testing_multiple_hybrids	Homo sapiens	Sushain Pandit	2009-01-02 00:00:00
testing_new_jar	Mus musculus	Sushain Pandit	2009-01-02 00:00:00
Gene expression profiles of fasting induced changes in liver and fat tissues of pigs expressing the MC4R D298N variant	Sus scrofa	Sender Lkhagvadorj	2008-09-19 00:00:00
Another Test by Neeraj	Homo sapiens	neeraj koul	2008-08-27 00:00:00
A test Experiment (Remove)	Homo sapiens	neeraj koul	2008-08-26 00:00:00
Expression data from non-infected and Salmonella Choleraesuis infected mesenteric lymph nodes	Sus scrofa	Keith Callenberg	2008-06-30 00:00:00
Analysis of differential gene expression in Salmonella shedder vs non-shedder pigs	Sus scrofa	Jolita Uthe	2008-04-02 00:00:00
Expression data from non-infected and Salmonella Typhimurium infected mesenteric lymph nodes	Sus scrofa	remy younes	2008-03-28 00:00:00

5.2 Regarding soft files and zip archives

Upon clicking on the experiment in Search Results, you will be presented by the experiment page with download links for soft and a zip archive.

- The soft file is accordance with GEO guidelines. This file refers other experiment-related files like .chp and .cel.
- In case of any missing valules that the system couldn't find, there would be a placeholder *-anexdb_admin insert value here*.

This field would indicate that there has been some missing information in the experiment submission and you should either change it to the correct value manually, or re-submit the experiment with the correct value / files.

- The zip archive would contain soft, .cel, .chp files for the experiment submission. All these would be needed since soft refer to them and they're all needed to conform to GEO submission requirements.
- It is highly recommended to cross-check the zip contents before submitting it to GEO. In case, any of the experiment-related files (.cel, .chp) is found missing, you should either manually add it to the archive, or check your experiment submission (on the experiments page) to ensure that you indeed submitted the file during the submission process.
- If you find that the file is absent in the zip even though it appears in the submission, please drop a note to one of the points of contact (refer 1.4).

5.3 Submitting the experiment to GEO

For successful submission to GEO, one would need to submit a zip archive of the generated SOFT file and the CEL/CHP files using the *Direct deposit* option at the [Direct Deposit](#) page.

- Select 'File format' SOFT, 'Submission kind' new and enter an appropriate 'Release date' for your data, before clicking the 'Submit' button.

5.4 Updating the soft after submission (if required)

For performing batch updates in the submitted SOFT, follow these guidelines: -

- Include the attribute "!Sample_GEO_accession = GSMxxx" where GSMxxx indicates the accession number of the record to be updated (similarly, use !Platform_GEO_accession = GPLxxx, and !Series_GEO_accession = GSExxx).
- You can provide the entire SOFT record with the necessary revisions. Alternatively, you can simply provide the revised attributes and/or data table. Data table updates are 'all-or-nothing', e.g., to revise data in one column of a table, it is necessary to update the whole table.
- Submit your SOFT update file by selecting the 'SOFT' option on the [Direct Deposit](#) page. Make sure to check the 'Update' box. Successful updates will be reflected immediately on your GEO records.

Refer the following examples (data tables, if any, truncated at 20 rows) for creating SOFT update files:

- An example of a file [updating Series attributes](#).
- An example of a file [updating an attribute in multiple Samples](#).
- An example of a file [updating data tables in multiple Samples](#).
- An example of a file [updating data tables and attributes in multiple Samples](#).
- An example of a file [updating Samples using complete SOFT records](#).